



ENDO RLCZ

Endoscopic Linear Cutting Staplers



**New
Design**



**Omnia Reload
Compatible**




**Improved
Force-to-fire**



CE Certificated




ENDO REACH RLCZ Universal Stapler

- Compatible with Omnia reloads
- Improved handle design



Product Code	Shaft Length (mm)	Max number of firing	QTY/BOX
ENDO RLCSZ	85	25	3
ENDO RLCZ	155		
ENDO RLCLZ	255		

Compatible Reload Specifications

Product Code	Staple Line Length (mm)	Trocar Compatibility (mm)	Color	Closed Staple Height (mm)	Tissue Thickness	Number of staple rows	QTY/BOX
ID30TAN	30	12		0.75/1.0/1.25	Vascular/ Thin	6	6
ID45TAN	45						
ID60TAN	60						
ID30PUL	30	15		1.25/1.5/1.75	Medium/ Thick		
ID45PUL	45						
ID60PUL	60						
ID45BLK	45	15		1.75/2.0/2.2	Extra Thick		
ID60BLK	60						

Contacts

Global Headquarter
120 Xinxing Road
West Zone, TEDA
300462 Tianjin
China

T +86-22-25326259
F +86-22-66200915

Importer



Distributor

