

REACH ACS-D Series

Circular Staplers



**New
Design**



**Staple height
control**



**Single-hand
Operation**



**Extended
shaft length**




CE Certificated

Enhanced Safety



Precise staple height control
with Staple Height Indicator window



More staples in the cartridge
providing better hemostasis with
maintained compression



Tactile and audio feedback
help physicians to confirm a
successful firing

Improved User Experience



Thinner anvil design
facilitates easy maneuvering in
tight spaces during procedures



180° arc Indicator window
easily confirm compression status
from multiple angles



Single-handle design
easier to maneuver and fire

Product Specification

Product Code	Color	Shaft Length	Outer Diameter	Lumen	Open Staple Height	Closed Staple Height	Number of Staples	Rows of Staples	QTY/BOX
ACS 21D	Gray	27 cm	21.5 mm	12.5 mm	3.5 mm	1.5-1.8 mm	16	2	3
ACS 25D	White		25.5 mm	16 mm	4.8 mm	1.5-2.2 mm	22		
ACS 28D	Blue		28.5 mm	18.8 mm			26		
ACS 31D	Green		31.5 mm	22.7 mm			30		
ACS 34D	Gold		34.5 mm	25 mm			32		

Contacts

Global Headquarter
120 Xinxing Road
West Zone, TEDA
300462 Tianjin
China

T +86-22-25326259
F +86-22-66200915

Imported by



Distributed by



This product is not available for Germany