



## Endo Reach ENDO III

Endoscopic Linear Cutting  
Staplers and Reloads



New Architecture



Improved  
compression



Single-hand  
operation



Single-use  
cartridge blade



CE Certified

## Endo Reach ENDO III Stapler

A compact design for maximum performance. Endo III stapling system is an ideal tool for surgical procedures, delivering consistent results with confidence.



**Improved compression**  
with new anvil design that maximizes compression and reduce tissue damage



**Smooth 45° articulating tip**  
45° articulating function in both directions



**Single-hand operation**  
Functional knob to easily switch between rotation and articulation.

Series	Product Code	Shaft Length	Max number of firings	QTY/BOX
<b>ENDO III 45mm Articulating Stapler</b>	REC45A S	282 mm	12	3
	REC45A	332 mm		
	REC45A L	432 mm		
<b>ENDO III 60mm Articulating Stapler</b>	REC60A S	298 mm		
	REC60A	348 mm		
	REC60A L	448 mm		
<b>ENDO III 45mm Articulating Stapler with Curved Tip</b>	REC45BA S	293 mm		
	REC45BA	343 mm		
	REC45BA L	443 mm		
<b>ENDO III 60mm Articulating Stapler with Curved Tip</b>	REC60BA S	309 mm		
	REC60BA	359 mm		
	REC60BA L	459 mm		

## Endo Reach ENDO III Cartridge

A wide range of Loading Units that offers excellent and smooth tissue dividing experience at 45 or 60 mm staple line length.



**Cartridge variety**  
6 staple heights in both 45 and 60 mm length are available for any surgical situation



**12mm-trocar accessibility**  
All cartridges, including green and black ones go through 12 mm trocars



**Single-use cartridge blade**  
Reduces the risk of cross contamination

Product Code	Staple Line Length	Cartridge Color	Open Staple Height	Closed Staple Height	Tissue Thickness	Number of staple rows	QTY/BOX
<b>REC 45 GRA</b>	45mm	Gray	2.0 mm	0.75 mm	Extra Thin	6	12
<b>REC 60 GRA</b>	60mm						
<b>REC 45 WHT</b>	45mm	White	2.5 mm	1.0 mm	Vascular		
<b>REC 60 WHT</b>	60mm						
<b>REC 45 BLU</b>	45mm	Blue	3.5 mm	1.5 mm	Thin		
<b>REC 60 BLU</b>	60mm						
<b>REC 45 GLD</b>	45mm	Gold	3.8 mm	1.75 mm	Medium		
<b>REC 60 GLD</b>	60mm						
<b>REC 45 GRN</b>	45mm	Green	4.1 mm	2.0 mm	Thick		
<b>REC 60 GRN</b>	60mm						
<b>REC 45 BLK</b>	45mm	Black	4.4 mm	2.2 mm	Extra Thick		
<b>REC 60 BLK</b>	60mm						

## Contacts

### Global Headquarter

120 Xinxing Road  
West Zone, TEDA  
300462 Tianjin  
China

**T** +86-22-25326259

**F** +86-22-66200915

## Importer



## Distributor

