



ENDO REACH

Endoscopic Linear Cutting Staplers



**New
Design**



**Omnia Reload
Compatible**




**Improved
Force-to-fire**

For LATAM Use Only
Available products subject to
local regulations and registrations

ENDO REACH Universal Stapler

- Compatible with Omnia reloads
- Improved handle design



Product Code	Shaft Length (mm)	Max number of firing	QTY/BOX
IDCS	85	25	3
IDCM	155		
IDCL	255		

Compatible Reload Specifications

Product Code	Staple Line Length (mm)	Trocar Compatibility (mm)	Color	Closed Staple Height (mm)	Tissue Thickness	Number of staple rows	QTY/BOX
ID30TAN IDSTV30TAN/IDSTV30TANB IDT30TAN-R/IDT30TAN-L	30	12	Tan	0.75/1.0/1.25	Vascular/ Thin	6	6
ID45TAN IDSTV45TAN/IDSTV45TANB IDT45TAN-R/IDT45TAN-L	45						
ID60TAN IDSTV60TAN/IDSTV60TANB IDT60TAN-R/IDT60TAN-L	60						
ID30PUL IDSTV30PUL/IDSTV30PULB IDT30PUL-R/IDT30PUL-L	30		Purple	1.25/1.5/1.75	Medium/ Thick		
ID45PUL IDSTV45PUL/IDSTV45PULB IDT45PUL-R/IDT45PUL-L	45						
ID60PUL IDSTV60PUL/IDSTV60PULB IDT60PUL-R/IDT60PUL-L	60						
ID45BLK/ IDSTV45BLK IDT45BLK-R/IDT45BLK-L	45	15	Black	1.75/2.0/2.2	Extra Thick		
ID60BLK/ IDSTV60BLK IDT60BLK-R/IDT60BLK-L	60						

Contacts

International Headquarters

16 Science Park Drive
#04-03 DNV GL
Technology Centre
Singapore 118227

T +65 6865 9879

Importer



Distributor

