



RLS

Linear Staplers with
Single Use Loading Units



30/45/60
Staple Length



Maximum
8 Firings



Superior Staple
Formation



Slim
Design

Product Features

- PST (Position Stapling Technology) provides reliable staple formation.
- Slim anvil improves access in narrow spaces.
- Parallel jaw closure provides consistent staple formation.
- Cutting edge provides ease of use.
- Slimmer and extended length improves access, especially in LAR(TME) procedures.
- Two-stage parallel closed design ensures easy tissue placement.

RLS Linear Stapler

Product Code	Stapler Line	Preloaded with unit	Max Firings	QTY/BOX
RLS3035S	30 mm	RLS3035L	8	3
RLS4535S	45 mm	RLS4535L		
RLS6035S	60 mm	RLS6035L		
RLS3048S	30 mm	RLS3048L		
RLS4548S	45 mm	RLS4548L		
RLS6048S	60 mm	RLS6048L		

RLS Linear Stapler Unit (SULU)

Product Code	Stapler Line	Cartridge Color	Open Staple Height	Closed Staple Height	Number of staples	Number of staple rows	QTY/BOX
RLS3035L	30 mm	Blue	3.5 mm	1.5 mm	11	2	6
RLS4535L	45 mm				15		
RLS6035L	60 mm				21		
RLS3048L	30 mm	Green	4.8 mm	2.0 mm	11		
RLS4548L	45 mm				15		
RLS6048L	60 mm				21		

202508

Global Headquarter

120 Xinxing Road
West Zone, TEDA
300462 Tianjin
China

T +86-22-25326259

F +86-22-66200915

